

## Mechanical Data

Maximum rotor speed	1760rpm
Maximum G-force	700g
Drive system	Direct drive
Maximum sample load	6 x 450ml / 5 liters

### Vacuum System

Pressure display	0-1200mbar
Pressure control	Automatic / 3mbar / to atmosphere
System ultimate vacuum	3mbar
Bumping / foaming protection	Dri-Pure*

### Temperature And Control

Control range	Ambient +7°C to 60°C
Control accuracy	±1°C
Temperature sensing	Via thermistor
Display range	0°C to 60°C
End of method	Time or automatic
Process visualization	Strobe & Delta T

### Solvent Compatibility

Boiling point range	40°C to 160°C at ambient
Includes	Alcohols, DCM/methylene chloride, DMF, ethyl acetate, water
HCl	Not compatible
Diethyl ether	Requires Inert Gas Purge option (compatible with flask rotor only)

### Dimensions

Width x Depth x Height	780 x 640 x 530mm
Height with lid open	782mm
Weight	75kg

### Services

Rocket Synergy 2 requires one mains power outlet	
UK & Europe	230V (±10%), 50Hz, 13A
USA	120V (±10%), 60Hz, 15A
Japan	100V (±10%), 50Hz or 60Hz, 15A
USB A	For data upload and download
For chamber water	Approx. 50ml per day

### Cold Trap Cooling Requirement

Temp range	-15°C to +10°C dependent upon application
Heat removal	700 Watts at +10°C
Flow rate	1.5 to 2.5 l/min
Pressure	1 (mini) to 2 bar (max) static
Connections (to chiller):	8mm nylon hardwall tube for Genevac supplied chiller

### Recirculating Chiller

A powerful recirculating chiller is available for the Rocket Synergy 2 evaporation system. The system can control the chiller via an RS232 link, thereby providing improved solvent recovery and better drying of samples compared with using a static cooled supply. A connection kit with insulated pipe work is available to accompany the chiller.



### Maintenance

All seals are durable consumables and user replaceable. Easy access is provided to the pump, which can be maintained by trained users.

### Safety

Complies with BS EN 61010-1:2010 and is CE marked.