

KEN BEDPAN WASHER DISINFECTOR 736

Intelligent Disinfection

- CLEAN, LEAN & GREEN

- √ Icons on display
- √ Total disinfecting system
- √ Easy to service access
- √ Complies with 15883

BWD 736 has been developed for hospitals, special wards and with its compact design it has a large capacity with a small depth.

The BWD 736 is easy to operate by touch panel with icons showing all relevant information of the

program needed by the operators.

The machine is functional with bedpans, urinal bottles, commodes, urine glasses or with options for cytostatica and clostridium difficile.

The different bedpan holders take all types of bedpans. The capacity of BWD 736 is three bedpans or 1 round bedpan and three urinal bottles per cycle. Most holders will be able to hold various commodes and urine measuring glasses.

A special utensils basket is available for washing kidney bowls and other small items.

All data from the BWD736 are fully exchangeable with common computer systems.





Advantages:

- Easy to use icons
- 7 standard programs
- Total disinfection guarantee
- Quick-start automatic closing of door and automatic programme start
- Door suitable as slopsink bedpan emptied automatically
 • Service access from the front / top

Hygiene safety

- EN 62366 tested
- EN ISO 15883-1-3-5

Environment

- Low energy costsLow noise level
- Low chemical usage

	VAP -
Technical description	Values
Outer dimensions / +plinth	H: 1540 x W: 600 x D: 460 mm
Operative chamber dimensions	H: 500 x W: 485 x D: 370 mm
Weight excl./incl. packaging	119 kg / + plinth 121 kg
Built in combitank	12 L
Built in detergent cabinet	2 x 5 L
Capacity bedpan / urinal bottles / bowls	3 pcs / 3 pcs / 3 pcs at a time
Bedpan program	7 min.
Sealing class	IP 21
Water connections	2 x 3/4 "
Cold water	200-800kPa, temp. 5-30°, 0-3° dH
Hot water	200-800kPa, temp. 5-30°, 0-3° dH
Sound pressure ISO 11204 L	<50 dB(A)
Water consumption per phase	See technical manual
Power supply options	3x200-415 V- 50/60 Hz
Heating elements in combitank	6 kW/ 9 kW
Programs available	7 standard
Drain	ø90 mm

Contact a KEN partner