

Portable Autoclave

Range of Application:

It's desk-top autoclave sterilizer which is designed for the hospital, clinic, lab, etc. It's mostly suitable for sterilizing of surgical, dental instruments, glassware, culture medium and biological dressing, food and goods, etc.

Features:

- Easy operation, safe and reliable.
- The sealing ring is made of silicon rubber.
- Two safety valves keep double safety (for timing control type).
- It has self-adjustment function (for timing control type).
- It is well sealed and effectively prevents steam leakage.
- Over temperature & over pressure auto protection.
- Electric heating or gas heating can be used (for Dual fuel Type).
- The sterilizer is provided with a power switch, user can control the operation easily.
- Immersion type electric heating pipe, thermal efficiency is high, stainless steel material tube is durable to use.
- Steam-water inner circulation system. No steam discharge and the environment for sterilizing will be clean and dry.
- Double groove pressure gauge, could clearly indicate, convenient observation, reliable performance.
- When at lower water level, or at the end of heat preservation time, have buzzer alarm, timing control type will automatically stop heating (for timing control type).



Timing control type



Anti-dry out type



Dual fuel type

Technical Parameters:

Model	BKM-P18(D)	BKM-P24(D)	BKM-P18(A)	BKM-P24(A)	BKM-P18(B)	BKM-P24(B)
Type	Timing control type		Anti-dry out type		Dual fuel Type	
Capacity	18L	24L	18L	24L	18L	24L
Chamber Size(mm)	φ280*210	φ280*310	φ280*210	φ280*310	φ280*210	φ280*380
Time Setting	1-99min		NO			
Temperature Setting	RT+5-126°C		126°C(default) or 121°C		126°C	
Consumption	20W				AC110/220v/19A, 60/50Hz	
Power Supply	AC220V, 50/60Hz				AC110/220v/19A, 60/50Hz	
Working Pressure	Along with temperature change		0.142Mpa		0.165Mpa	
Package Size(W*D*H)mm	420*420*650	420*420*740	420*420*650	420*420*740	420*420*470	410*410*550
Gross Weight(kg)	12	14	12	14	8	17